



# Pleural Fluid (Biochemistry)

## Synonyms

## Clinical Indication

Pleural effusion. To determine whether the fluid is a transudate or exudate.  
Serum sample taken within 24 hours of the fluid sample is required.  
Refer to Fluid (biochemistry) for ?chylothorax/pseudochylothorax information.

## Part of Profile / See Also

Fluid total protein, cholesterol, LDH and triglycerides are measured.

## Request Form

Combined Pathology manual Blood form or ICE request

## Availability / Frequency of Analysis

On request

## Turnaround Time

1 day

## Patient Preparation

None required

## Sample Requirements

### Specimen Type


Fluid. Serum sample drawn 24 hours pre or post fluid collection

### Volume

1ml of fluid in White top universal.


### Required Containers

**This should contain Fluid.**



White top universal container (total protein, cholesterol, triglycerides and LDH.)

**This should contain Blood**



Yellow top (SST) tube

## Reference Range & Units

A fluid score is calculated from fluid total protein, cholesterol and LDH levels. A score of <0.75 is consistent with a transudate; >1.5 is consistent with an exudate. (Further details on in-house study available)

## Interferences

Haemolysis.

## Interpretation & Clinical

### Decision Value (if applicable)

For investigation of suspected chylothorax, a fluid with a triglyceride level above 1.24 mmol/L has a 99% chance of being chylous whereas a fluid with a triglyceride level below 0.57 mmol/L has no more than 5% chance.

## References

**Test code**

PLEU (Fluid LDH and Fluid Protein). FCHO (Fluid cholesterol). FTG (Fluid triglyceride)

**Lab Handling**

Record appearance before centrifugation.  
Centrifuge universal and store at -20°C and store aliquot at 4°C.