



Accredited to
ISO 15189:2022

Triglycerides

Synonyms

Clinical Indication

High serum triglyceride concentrations are probably an independent (of cholesterol) risk factor for coronary heart disease. Very high levels (>10 mmol/L) confer a risk of pancreatitis.

Triglycerides circulate in blood as components of various lipoproteins, particularly very light density lipoproteins (VLDL) and chylomicrons.

If drug treatment is considered appropriate in hypercholesterolaemia, the choice of drug may be influenced by the additional presence of high triglyceride levels

Part of Profile / See Also

Lipids (fasting lipid profile)

Request Form

Combined Pathology manual Blood form or ICE request

Availability / Frequency of Analysis

On request.

Turnaround Time

Same day

Patient Preparation

Patients must fast for 10 hours prior to blood collection (for follow up tests diabetic patients may not need to fast, but levels may be increased as a result of recent food intake).

Please note: Samples are still analysed if the patient is not fasting.

Sample Requirements

Specimen Type

Serum and plasma

Volume

2 ml

Acceptable Containers



Yellow top (SST) tube



Green top (lithium-heparin) tube



paediatric orange top (lithium-heparin)



paediatric green top (lithium-heparin)

Plain serum samples may also be used.

Reference Range & Units

Fasting: Less than 2.26 mmol/L

Reference: National Cholesterol Education Program (NCEP) Unit Conversion

Interferences
Interpretation & Clinical
Decision Value (if applicable)

High triglyceride levels may be genetic or acquired due to drugs e.g. thiazides, corticosteroids, diabetes, obesity, or excess alcohol. Elevated triglyceride levels puts patients at increased risk of pancreatitis. Amylase may be added to the sample by the duty biochemist at Clinical authorisation

Critical Difference: 0.9 mmol/L

References

Beckman kit insert

Test code

TRIG LIP (fasting lipid profile)

Lab Handling

Analysed from primary tube and stored at 4°C.

Serum and plasma stable for 7 days at 4°C.