

## Farmers' Lung

Synonyms	Hypersensitivity pneumonitis (HP) or extrinsic allergic alveolitis.
Clinical Indication	A serological test to Thermophilic actinomycetes to aid diagnosis of Farmers Lung.
	Please provide detailed clinical information of occupational risk.
Part of Profile / See Also	
Request Form	ICE order comms form (Handwritten request – only if order comms not available).
Availability / Frequency of Analysis	Referred test: Available on request during normal working hours. Sent to Bristol Mycology Reference Laboratory (8043) overnight by courier.
Turnaround Time	10 days.
Patient Preparation	None
Sample Requirements	
Specimen Type	Blood
Volume	5ml
Container	
	Yellow top (SST) tube
	MICROTAINER
	Paediatric Yellow top (SST) tube
	Samples should be transported to laboratory immediately.
Reference Range & Units	
Interferences	Grossly haemolysed, contaminated or highly lipaemic samples.
Interpretation & Clinical	The reference laboratory provides interpretive comments which will be
Decision Value (if applicable)	included in the final report.
References	
Test code	FUNS
Lab Handling	Aliquot and store at 4°C; Send sample to reference laboratory as soon as possible.



PF-PTD-134



Accredited to ISO 15189:2022