



HFE Genotyping

Synonyms	Haemochromatosis genetics
Clinical Indication	For confirmation of diagnosis of hereditary haemochromatosis when transferrin saturation >50% or requested by Consultant Haematologist.
Part of Profile / See Also	
Request Form	Combined Pathology manual blood form or ICE request
Availability / Frequency of Analysis	Referred test: Analysed by Barts Health if specific criteria met (8285)
Turnaround Time	7 - 10 days
Patient Preparation	None required
Sample Requirements	
Specimen Type	Whole blood (EDTA)
Volume	1 x 3 or 4 mL
Container	 Purple top (EDTA) tube.
	 Paediatric Lavender top (EDTA) tube
Reference Range & Units	N/A
Interferences	
Interpretation & Clinical Decision Value (if applicable)	See individual reports.
References	
Test code	HFE
Lab Handling	Samples should be sent to referral laboratory by post or courier.