

Culture for aerobic and anaerobic bacteria.

handwritten Microbiology request form



PF-PTD-238

Pus

Synonyms

Clinical Indication

Part of Profile / See Also

Request Form

Availability / Frequency of

Analysis

. . .

During working day

For urgent requests please contact laboratory on 01268 968285

Refrigerate sample if there is to be a delay in delivery to laboratory.

indicated or specifically requested. Refer to PF-PTD 215

Please give detailed clinical information, including date of onset/contact with the suspected infection, and any current, or intended, antibiotic therapy

Microscopy for acid fast bacilli (AFB) and culture for Mycobacteria if clinically

ICE / Medway order comms request or if order coms not available, use

Microscopy: For organisms and WBCs (parasites only where indicated)

Turnaround Time

Patient Preparation

Sample Requirements

Specimen Type

Fluid pus is preferred; a bacterial swab in transport medium is a poor

alternative.

10 days

Volume Up to 20 ml (minimum volume 0.5 ml)

Container White capped sterile universal - CE marked

Reference Range & Units

Interferences

Interpretation & Clinical

Decision Value (if applicable)

Cultures will be reported with antibiotic susceptibilities where appropriate.

ReferencesUK SMI Bacteriology | B 14 | Issue no: 6.2 | Issue date: 22.11.16 | Investigation

of Pus and Exudates

Test code PUSF

Lab Handling Process on receipt





PF-PTD-238