





8285
Accredited to
ISO 15189:2012

PNH

Synonyms	Paroxysmal Nocturnal Haemoglobinuria
Clinical Indication	PNH suspected by Haematology consultant. Clinically indicated if requested by Consultant Haematologist or if agreed with Consultant Haematologist.
Part of Profile / See Also	CD14 and CD16
Request Form	Combined Pathology manual blood form or ICE request
Availability / Frequency of Analysis	Monday to Thursday. Urgent samples only on Friday. Referred test: Analysed by Barts Health
Turnaround Time	7 days
Patient Preparation	None required
Sample Requirements	
Specimen Type	Whole blood
Volume	1 EDTA sample
Container	 Purple top (EDTA) tube.  Paediatric Lavender top (EDTA) tube
	Sample must be transported to the laboratory immediately.
Reference Range & Units	See individual report
Interferences	Prolonged storage
Interpretation & Clinical Decision Value (if applicable)	Immunophenotyping results are interpreted alongside other haematology results including FBC and film.
References	Barts Health SIHMDS User guide
Test code	PNH
Lab Handling	Send by courier, preferably by Thursday run, to the referral laboratory.