



PF-PTD-315

Wound Swab

Synonyms

Clinical Indication

If fluid pus is available, for example, from an abscess, a sample of this should be submitted in preference to a swab, see pus in test directory for more information

Please give detailed clinical information, including date of onset/contact with the suspected infection and any current or intended antibiotic therapy

Part of Profile / See Also

Microscopy: dependent on site and clinical information received.

Culture: all wound swabs are cultured for aerobic organisms. Culture for anaerobic organisms and other investigations may be performed, dependent on sample site and clinical information received.

Request Form

Order comms request i.e., ICE or Medway

or clearly handwritten Microbiology forms where order comms are not available

Availability / Frequency of

Analysis

On request during normal working hours

Bacterial swab in transport medium.

Turnaround Time

5 days

Patient Preparation

Sample Requirements

Samples should be transported to the laboratory on the day of collection.

Specimen Type

Volume

Container

As above

Reference Range & Units

Interferences

Quality of the sample may deteriorate if there is a delay in processing and isolation of certain organisms particularly anaerobes maybe compromised.

Interpretation & Clinical

Decision Value (if applicable)

References

UK SMI Bacteriology | B 11 | Issue no: 6.5 | Issue date: 19.12.18 |

Investigation of skin and superficial soft tissue infections

Test code

RCUL

Lab Handling

Samples should be processed as soon as possible.