

LFT, Bone

On request.

Same day

Not measured directly.



PF-PTD-354

Globulin

Accredited to ISO 15189:2012

Synonyms

Clinical Indication Investigation of protein abnormalities

Part of Profile / See Also

Request Form

Availability / Frequency of

Analysis

Turnaround Time

Patient Preparation

Sample Requirements

Specimen Type

Volume

Container

Reference Range & Units

20 - 35 g/L

If globulin level is raised and there is no history of previous investigation, serum electrophoresis will automatically be performed by the laboratory

Combined Pathology manual blood request form or ICE request

Calculation: Total protein-Albumin = Globulin g/L

Interferences

Interpretation & Clinical

Decision Value (if applicable)

References

Test code

L (part of LFT Profile) or BONE (part of bone profile)

Lab Handling

Analysed from the primary tube and stored at 4°C