



PF-PTD-54

Blood Film

Synonyms

Blood film, blood smear

Clinical Indication

Checking of morphological features of blood cells. Not to be used for malarial film request, please see Malarial Parasites.

Must always accompany a Full Blood Count. Reason for request should be given.

Part of Profile / See Also

FBC

Request Form

Combined Pathology manual Blood form or ICE request

Availability / Frequency of

On request / generated by Haematology as part of a full blood count if required.

Analysis

Routine request: 48 hours from time sample received in the laboratory.

Patient Preparation

Turnaround Time

None required

Sample Requirements

Specimen Type Whole Blood

Volume 3 or 4 mL

Container



Purple top (EDTA) tube



Paediatric Lavender top (ETDA) tube

Sample taken for FBC can also be used for blood film review. Sample should ideally reach the laboratory within 4 hours of collection to minimise effects of storage artefacts.

Reference Range & Units

Not applicable

Interferences

Clotted samples unsuitable, EDTA changes due to morphological features (see

sample requirements)

Interpretation & Clinical

Decision Value (if applicable)

Film morphology comments reported by laboratory, to be interpreted by clinicians.

Films will be passed onto the haematology clinicians for clinical comments on request of if clinical referral criteria are met.

References

N/A

Test code

FILM

Lab Handling

Prepare film as soon as possible.